

Detailed Hydrocarbon Analysis (DHA)

-Paraffins mix

1. Analytical Condition

YL 6500 Series

Oven : 35°C(15min)->1°C/min->60°C(20min)->2°C/min->200°C(5min)

Column : HP-5MS(15m*0.25mm*1.0um)/CP-Sil PONA CB(100m*0.25mm*0.5um)

Carrier gas : He, 2ml/min (SR 150:1)

Injector : Cap 250°C

Detector : FID 270°C (Mkup 20ml/min)

Injection Volume : 0.2ul

2. Chromatogram

