

Detailed Hydrocarbon Analysis (DHA)

-Paraffins mix

1. Analytical Condition

YL 6500 Series

Oven : Cryogenic, 5°C(10min)->5°C/min->50°C(0.1min)->1.5°C/min->200°C(10min)

Column : RTX-100-DHA (100m*0.25mm*0.5um)

Carrier gas : He, 2ml/min (SR 150:1)

Injector : Cap 250°C

Detector : FID 250°C (Mkup 20ml/min)

Injection Volume : 0.2ul

2. Chromatogram

