



## Diethyl phthalate & Transfluthrin Analysis

### 1. Analytical Condition

#### **YL 6500 Series**

Oven : 100°C (for 1 min)- temp-programed 10 °C/min-270 °C (for 1 min)

Column : HP-1 ( 30m\*0.32mm\*1.0um )

Carrier gas : He, 3ml/min (Split ratio 10:1)

Injector : Capillary 270 °C

Detector : FID 270 °C

Injection Volume : 1.0ul (Liquid)

### 2. Chromatogram

