

## Ethanol in Hand Cleanser

### 1. Analytical Condition

#### YL6500 GC

Oven : 50°C ( 1 min ) -> 20 °C/min -> 150 °C (0.1min)

Column : HP-INNOWAX ( 30 m \* 0.53 mm \* 1 µm)

Carrier gas : He, 8mL/min (Split ratio 20:1)

Injector : Capillary 230 °C

Detector : FID 250 °C

Injection Volume: 1 µL (Diluted in DMF)

### 2. Chromatogram

