

CCl(F)Cl

## HCFC-141b Purity Analysis (1,1-Dichloro-1-fluoroethane)

### 1. Analytical Condition

#### YL 6500 Series

Oven : 50°C (for 5 min)- temp-programed 15 °C/min-150 °C (for 1 min)

Column : GC-IC Wax (30m\*0.53mm\*1.00um)

Carrier gas : N2, 6ml/min (Split ratio 30:1)

Injector : Capillary 230°C

Detector : FID 230°C

Injection Volume : 1.0ul (Liquid)

### 2. Chromatogram

