

Acrylate Analysis

1. Analytical Condition

YL 9100 HPLC

Column : Polaris dC18 (4.6*250mm, 5um)

Mobile Phase : A(Water, 0.1% H3PO4)

B(ACN, 0.1% H3PO4)

Detector : UV/Vis 220nm

Flow : 1.0ml/min

Injection volume : 15 μ l

Oven Temp : 40°C

< Gradient Program >

Time(min)	A	B
init	100	0
1	100	0
11	0	100
16	0	100
16.1	100	0

2. Chromatogram

