



1,2-Benzisothiazolinone-3-one Analysis

1. Analytical Condition

YL9100 HPLC

Mobile Phase : Buffer(1) : ACN = 75 : 25

*Buffer(1) : [10mM NaH₂PO₄+10mM Na₂HPO₄+1-Octane sulfonic acid Na 700mg
+ 1ml phosphoric acid]/ 1L Water

Flow Rate : 1ml/min

Column Oven : 40°C

Detector : UV/Vis 254nm

Column : C18 (4.6*250mm, 5um)

Injection Volume : 20ul

2. Chromatogram

