

Citric acid Analysis

In Plating solutions of Electronic Component

1. Analytical Condition

YL 9100 Series

Mobile phase: ACN / 0.1% H₃PO₄ in water

Flow rate : 1 mL/min

Oven : 40°C

Detector : UV 210 nm

Column : SB-Aq C18(4.6 * 250 mm, 5 um)

Injection : 20 uL

Time (min.)	ACN	0.1% H ₃ PO ₄
0	0	100
7.00	0	100
7.10	100	0
17.00	100	0
17.1	0	100

2. Chromatogram

