

Ginkgo Biloba Leaf Extract Analysis - Ginkgolide

1. Analytical Condition

YL9100 HPLC

Mobile Phase : Water, MeOH

Flow Rate : 1.0ml/min

Column Oven : 40°C


Detector : ELSD 60°C, 2L/min

Column : C18 (4.6*250mm, 5um)

Injection Volume : 10ul

Gradient Program	Water	MeOH
Initial	75	25
23.00	52	48
30.00	10	90
33.00	10	90
34.00	75	25
45.00	75	25

2. Chromatogram



1. Bilobalide
2. Ginkgolide C
3. Ginkgolide A
4. Ginkgolide B

