

Salicylic acid, Pyridoxin HCL,

Resorcinol,

Mononitroguajacol Sodium Analysis

1. Analytical Condition

YL9100 HPLC

Mobile Phase : aqueous solution 700ml (10mM NaH₂PO₄+10mM Na₂HPO₄) + MeOH 300ml
+ 1g 1-Octane sulfonic acid Na => Mix at 1 bottle

Flow Rate : 1.0 ml/min

Column Oven : 35°C

Detector : UV/Vis 280 nm

Column : C18 (4.6*150mm, 5um)

Injection Volume : 10ul

2. Chromatogram

