

Parabens 5mix Analysis

1. Analytical Condition

	Time(min)	Water(0.1%Acetic acid)	ACN
YL 9100 HPLC Column : SP C18 MGII (4.6*250mm, 5um) Detector : PDA 257nm Flow : 1.0ml/min Injection volume : 20 μ l Oven Temp : 35	init	55	45
	5	50	50
	10	45	55
	15	40	60
	20	40	60
	23	35	65
	28	35	65
	30	0	100
	36	0	100
	37	55	45

2. Chromatogram

