

Parabens Analysis

1. Analytical Condition

YL 9100 HPLC

Column : SP C18 MG II
(4.6*250mm, 5um)

Detector : UV/Vis or PDA 257nm

Flow : 1.0ml/min

Injection volume : 20 μ l

Oven Temp : 35°C

Time(min)	Water(0.1%Acetic acid)	ACN
init	55	45
5	50	50
10	45	55
15	40	60
20	40	60
23	35	65
28	35	65
30	0	100
31	55	45

2. Chromatogram

