

Flavonoid Analysis

1. Analytical Condition

YL9100 HPLC

Mobile Phase : A = Water(0.1% H₃PO₄),

B = ACN(0.1% H₃PO₄)

Flow Rate : 0.8ml/min

Column Oven : 40°C

Detector : UV/Vis 280nm

Column : C18 (4.6*150mm, 5um)

Injection Volume : 20ul

<Gradient Program>

	A	B
0min	80	20
2min	80	20
55min	45	55
56min	80	20

2. Chromatogram

