

PMMA(Poly(Methyl Methacrylate)) Standard

1. Analytical Condition

YL GPC Analyzer

Mobile phase : 50mM LiBr in DMF

Flow rate : 1 mL / min

Injection volume : 100 $\mu\ell$

Column Temp : 35 $^{\circ}\text{C}$

Detector : RID, 35 $^{\circ}\text{C}$

Column : Waters Styragel HR 1 / 4E / 5E (7.8 * 300 mm)

+ Guard column

Standard

STD : PMMA, MW 54,500 / 31,600 / 10,100 / 2,580

2. Chromatogram

