

# Carbamate Pesticides(11 compounds) in Agricultural Products

## 1. Analytical Condition

### YL Carbamate Analyzer -FLD (OPA Derivatization)

Column : Carbamate C18 (4.6 X 250mm,5 $\mu$ m)

Mobile Phase : Water : AC N= Gradient

Detector : FLD Ex 340nm/Em 455nm

Injection : 20 $\mu$ l

Column Oven : 40°C

<Pinnacle PCX with 0.5ml reactor>

## 2. Chromatogram

